


**Test Sieves, 203 x 50 mm****Retsch****Test Sieves, 203 x 50 mm****Features:**

- Previously unattained product quality owing to advanced manufacturing processes and continuous optical inspection.
 - High quality joining technology guarantees permanently taut wire mesh.
 - High level of corrosion resistance and easy to clean owing to the use of high-alloy stainless steel.
 - "one-piece cast" sieve frame and seamless mesh joint prevent cross contamination.
 - Maximum stability and optimum density of the sieve stack.
 - 15% lighter than traditional sieves whilst at the same time having larger sieving area.
 - Clear sieve marking and traceability by individual laser engraving on each sieve.
- DIN ISO 3310-1. Dimensions (dia. x H): 203mm x 50mm. Other sieve sizes available on request



Code	Description	Packaging
LLG09226865	Test Sieves, 203 x 50 mm, Mesh size 100 µm	1 pz.
LLG09226866	Test Sieves, 203 x 50 mm, Mesh size 106 µm	1 pz.
LLG09226867	Test Sieves, 203 x 50 mm, Mesh size 112 µm	1 pz.
LLG09226868	Test Sieves, 203 x 50 mm, Mesh size 125 µm	1 pz.
LLG09226869	Test Sieves, 203 x 50 mm, Mesh size 140 µm	1 pz.
LLG09226870	Test Sieves, 203 x 50 mm, Mesh size 150 µm	1 pz.
LLG09226871	Test Sieves, 203 x 50 mm, Mesh size 160 µm	1 pz.
LLG09226872	Test Sieves, 203 x 50 mm, Mesh size 180 µm	1 pz.
LLG09226850	Test Sieves, 203 x 50 mm, Mesh size 20 µm	1 pz.
LLG09226873	Test Sieves, 203 x 50 mm, Mesh size 200 µm	1 pz.
LLG09226874	Test Sieves, 203 x 50 mm, Mesh size 212 µm	1 pz.
LLG09226875	Test Sieves, 203 x 50 mm, Mesh size 224 µm	1 pz.
LLG09226851	Test Sieves, 203 x 50 mm, Mesh size 25 µm	1 pz.
LLG09226876	Test Sieves, 203 x 50 mm, Mesh size 250 µm	1 pz.
LLG09226877	Test Sieves, 203 x 50 mm, Mesh size 280 µm	1 pz.
LLG09226878	Test Sieves, 203 x 50 mm, Mesh size 300 µm	1 pz.
LLG09226879	Test Sieves, 203 x 50 mm, Mesh size 315 µm	1 pz.
LLG09226852	Test Sieves, 203 x 50 mm, Mesh size 32 µm	1 pz.
LLG09226880	Test Sieves, 203 x 50 mm, Mesh size 355 µm	1 pz.
LLG09226853	Test Sieves, 203 x 50 mm, Mesh size 36 µm	1 pz.



Code	Description	Packaging
LLG09226854	Test Sieves, 203 x 50 mm, Mesh size 38 µm	1 pz.
LLG09226855	Test Sieves, 203 x 50 mm, Mesh size 40 µm	1 pz.
LLG09226881	Test Sieves, 203 x 50 mm, Mesh size 400 µm	1 pz.
LLG09226882	Test Sieves, 203 x 50 mm, Mesh size 425 µm	1 pz.
LLG09226856	Test Sieves, 203 x 50 mm, Mesh size 45 µm	1 pz.
LLG09226883	Test Sieves, 203 x 50 mm, Mesh size 450 µm	1 pz.
LLG09226857	Test Sieves, 203 x 50 mm, Mesh size 50 µm	1 pz.
LLG09226884	Test Sieves, 203 x 50 mm, Mesh size 500 µm	1 pz.
LLG09226858	Test Sieves, 203 x 50 mm, Mesh size 53 µm	1 pz.
LLG09226859	Test Sieves, 203 x 50 mm, Mesh size 56 µm	1 pz.
LLG09226885	Test Sieves, 203 x 50 mm, Mesh size 560 µm	1 pz.
LLG09226886	Test Sieves, 203 x 50 mm, Mesh size 600 µm	1 pz.
LLG09226860	Test Sieves, 203 x 50 mm, Mesh size 63 µm	1 pz.
LLG09226887	Test Sieves, 203 x 50 mm, Mesh size 630 µm	1 pz.
LLG09226861	Test Sieves, 203 x 50 mm, Mesh size 71 µm	1 pz.
LLG09226888	Test Sieves, 203 x 50 mm, Mesh size 710 µm	1 pz.
LLG09226862	Test Sieves, 203 x 50 mm, Mesh size 75 µm	1 pz.
LLG09226863	Test Sieves, 203 x 50 mm, Mesh size 80 µm	1 pz.
LLG09226889	Test Sieves, 203 x 50 mm, Mesh size 800 µm	1 pz.
LLG09226890	Test Sieves, 203 x 50 mm, Mesh size 850 µm	1 pz.



Code	Description	Packaging
LLG09226864	Test Sieves, 203 x 50 mm, Mesh size 90 µm	1 pz.
LLG09226891	Test Sieves, 203 x 50 mm, Mesh size 900 µm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

