

**Racks for cryoboxes Arctic Squares, Stainless steel****Crystal Technology &****Racks for cryoboxes Arctic Squares®, Stainless steel**

Stainless steel. Fits most standard upright freezers. Overall rack size (H x W x D): 238 x 139 x 571 mm. With four drawers to hold 16 boxes of 133 x 133 x 53 mm. Drawers feature a unique stop that when fully extended pivot down allowing for easy box access. Freezer rack also includes a pull-out wire handle for easy removal off the shelf. Polycarbonate drawer handles provide a user-friendly grip. Includes a slot and label for drawer/sample identification.



| Code | Description | Packaging |
|-------------|--|-----------|
| LLG09193980 | Racks for cryoboxes Arctic Squares®, Stainless steel, Type Arctic Squares® | 1 pz. |



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

