



Conductivity Standards

Conductivity Standards

Conductivity standards tested and ISO 17025 certified at 25 °C. Uncertainty of measurement of $\pm 1\%$. Directly traceable to N.I.S.T. Temperature dependence data is printed on the label as are lot numbers and expiry dates. Supplied in 500 ml HDPE bottle. 1.3 $\mu\text{S}/\text{cm}$ in 250 ml bottle. Calibration values 84, 147, 1413, 12880 $\mu\text{S}/\text{cm}$. Control Values 1.3 - 350000 $\mu\text{S}/\text{cm}$.



Code	Description	Packaging
LLG06225059	Conductivity Standards, Type 1.3* $\mu\text{S/cm}$, 250 ml	1 pz.
LLG06074754	Conductivity Standards, Type 50 $\mu\text{S/cm}$, 500 ml	1 pz.
LLG06205412	Conductivity Standards, Type 5000 $\mu\text{S/cm}$, 500 ml	1 pz.
LLG06205413	Conductivity Standards, Type 100000 $\mu\text{S/cm}$, 500 ml	1 pz.
LLG06205415	Conductivity Standards, Type 500 $\mu\text{S/cm}$, 500 ml	1 pz.
LLG06210129	Conductivity Standards, Type 50000 $\mu\text{S/cm}$, 500 ml	1 pz.
LLG06234859	Conductivity Standards, Type 100 $\mu\text{S/cm}$, 500 ml	1 pz.
LLG06802901	Conductivity Standards, Type 20 $\mu\text{S/cm}$, 500 ml	1 pz.
LLG06802902	Conductivity Standards, Type 200 $\mu\text{S/cm}$, 500 ml	1 pz.
LLG06901622	Conductivity Standards, Type 10 $\mu\text{S/cm}$, 500 ml	1 pz.
LLG06901670	Conductivity Standards, Type 1000 $\mu\text{S/cm}$, 500 ml	1 pz.
LLG06901671	Conductivity Standards, Type 10000 $\mu\text{S/cm}$, 500 ml	1 pz.
LLG09040904	Conductivity Standards, Type 84 $\mu\text{S/cm}$, 500 ml	1 pz.
LLG09040905	Conductivity Standards, Type 147 $\mu\text{S/cm}$, 500 ml	1 pz.
LLG09040906	Conductivity Standards, Type 1413 $\mu\text{S/cm}$, 500 ml	1 pz.
LLG09040907	Conductivity Standards, Type 12880 $\mu\text{S/cm}$, 500 ml	1 pz.
LLG09040987	Conductivity Standards, Type 350000 $\mu\text{S/cm}$, 500 ml	1 pz.
LLG06901629	Conductivity Standards, Type 5 $\mu\text{S/cm}$, 500 ml	1 pz.





CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

