



## Conductivity Standards

### Conductivity Standards

Conductivity standards tested and ISO 17025 certified at 25 °C. Uncertainty of measurement of  $\pm 1\%$ . Directly traceable to N.I.S.T. Temperature dependence data is printed on the label as are lot numbers and expiry dates. Supplied in 500 ml HDPE bottle. 1.3  $\mu\text{S}/\text{cm}$  in 250 ml bottle. Calibration values 84, 147, 1413, 12880  $\mu\text{S}/\text{cm}$ . Control Values 1.3 - 350000  $\mu\text{S}/\text{cm}$ .



| Code        | Description                                                   | Packaging |
|-------------|---------------------------------------------------------------|-----------|
| LLG06225059 | Conductivity Standards, Type 1.3* $\mu\text{S/cm}$ , 250 ml   | 1 pz.     |
| LLG06074754 | Conductivity Standards, Type 50 $\mu\text{S/cm}$ , 500 ml     | 1 pz.     |
| LLG06205412 | Conductivity Standards, Type 5000 $\mu\text{S/cm}$ , 500 ml   | 1 pz.     |
| LLG06205413 | Conductivity Standards, Type 100000 $\mu\text{S/cm}$ , 500 ml | 1 pz.     |
| LLG06205415 | Conductivity Standards, Type 500 $\mu\text{S/cm}$ , 500 ml    | 1 pz.     |
| LLG06210129 | Conductivity Standards, Type 50000 $\mu\text{S/cm}$ , 500 ml  | 1 pz.     |
| LLG06234859 | Conductivity Standards, Type 100 $\mu\text{S/cm}$ , 500 ml    | 1 pz.     |
| LLG06802901 | Conductivity Standards, Type 20 $\mu\text{S/cm}$ , 500 ml     | 1 pz.     |
| LLG06802902 | Conductivity Standards, Type 200 $\mu\text{S/cm}$ , 500 ml    | 1 pz.     |
| LLG06901622 | Conductivity Standards, Type 10 $\mu\text{S/cm}$ , 500 ml     | 1 pz.     |
| LLG06901670 | Conductivity Standards, Type 1000 $\mu\text{S/cm}$ , 500 ml   | 1 pz.     |
| LLG06901671 | Conductivity Standards, Type 10000 $\mu\text{S/cm}$ , 500 ml  | 1 pz.     |
| LLG09040904 | Conductivity Standards, Type 84 $\mu\text{S/cm}$ , 500 ml     | 1 pz.     |
| LLG09040905 | Conductivity Standards, Type 147 $\mu\text{S/cm}$ , 500 ml    | 1 pz.     |
| LLG09040906 | Conductivity Standards, Type 1413 $\mu\text{S/cm}$ , 500 ml   | 1 pz.     |
| LLG09040907 | Conductivity Standards, Type 12880 $\mu\text{S/cm}$ , 500 ml  | 1 pz.     |
| LLG09040987 | Conductivity Standards, Type 350000 $\mu\text{S/cm}$ , 500 ml | 1 pz.     |
| LLG06901629 | Conductivity Standards, Type 5 $\mu\text{S/cm}$ , 500 ml      | 1 pz.     |





CARLO ERBA Reagents S.A.S.  
Chaussée du Vexin – 27106 Val  
de Reuil cedex  
N° TVA: FR 63391048824  
Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



*CARLO ERBA  
Reagents operates with  
a Certified Quality  
Management System*

