



Cleaning and storage solutions for pH-electrodes

Cleaning and storage solutions for pH-electrodes

Cleaning, rinse, regeneration and storage solution. Designed to extend the useful life of pH electrodes.



Code	Description	Packaging
LLG06287444	Cleaning and storage solutions for pH-electrodes, Type Electrode Rinse solution, 500 ml	1 pz
LLG09040940	Cleaning and storage solutions for pH-electrodes, Type Pepsin / HCl, for protein removal, 500 ml	1 pz
LLG09040941	Cleaning and storage solutions for pH-electrodes, Type Thiourea / HCl, for sulphide removal, 500 ml	1 pz
LLG09040942	Cleaning and storage solutions for pH-electrodes, Type Organic cleaning solution, 500 ml	1 pz
LLG09040943	Cleaning and storage solutions for pH-electrodes, Type Electrode Regeneration solution, 100 ml	1 pz
LLG09040944	Cleaning and storage solutions for pH-electrodes, Type Electrode Storage solution, 500 ml	1 pz
LLG09040945	Cleaning and storage solutions for pH-electrodes, Type Electrode Storage solution, 100 ml	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

