



## Cleaning and storage solutions for pH-electrodes

### Cleaning and storage solutions for pH-electrodes

Cleaning, rinse, regeneration and storage solution. Designed to extend the useful life of pH electrodes.



Code	Description	Packaging
LLG09040943	Cleaning and storage solutions for pH-electrodes, Type Electrode Regeneration solution, 100 ml	1 pz.
LLG09040945	Cleaning and storage solutions for pH-electrodes, Type Electrode Storage solution, 100 ml	1 pz.
LLG09040940	Cleaning and storage solutions for pH-electrodes, Type Pepsin / HCl, for protein removal, 500 ml	1 pz.
LLG09040941	Cleaning and storage solutions for pH-electrodes, Type Thiourea / HCl, for sulphide removal, 500 ml	1 pz.
LLG09040942	Cleaning and storage solutions for pH-electrodes, Type Organic cleaning solution, 500 ml	1 pz.
LLG09040944	Cleaning and storage solutions for pH-electrodes, Type Electrode Storage solution, 500 ml	1 pz.
LLG06287444	Cleaning and storage solutions for pH-electrodes, Type Electrode Rinse solution, 500 ml	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

