



Metal caps, Labocap without handle

schuett-biotec

Metal caps, Labocap without handle



Code	Description	Packaging
LLG09231423	Metal caps, Labocap without handle, For neck diam. 14/15 mm, Colour blue, Handle without	100 pz.
LLG09231385	Metal caps, Labocap without handle, For neck diam. 24/26 mm, Colour silver, Handle without	100 pz.
LLG09231428	Metal caps, Labocap without handle, For neck diam. 24/26 mm, Colour blue, Handle without	100 pz.
LLG09231448	Metal caps, Labocap without handle, For neck diam. 24/26 mm, Colour red, Handle without	100 pz.
LLG09231449	Metal caps, Labocap without handle, For neck diam. 28/30 mm, Colour red, Handle without	100 pz.
LLG09231334	Metal caps, Labocap without handle, For neck diam. 12/13 mm, Colour black, Handle without	100 pz.
LLG09231337	Metal caps, Labocap without handle, For neck diam. 17/18 mm, Colour black, Handle without	100 pz.
LLG09231341	Metal caps, Labocap without handle, For neck diam. 19/20 mm, Colour black, Handle without	100 pz.
LLG09231466	Metal caps, Labocap without handle, For neck diam. 15/16 mm, Colour black, Handle without	100 pz.
LLG09231421	Metal caps, Labocap without handle, For neck diam. 9/10 mm, Colour blue, Handle without	100 pz.
LLG09231422	Metal caps, Labocap without handle, For neck diam. 12/13 mm, Colour blue, Handle without	100 pz.
LLG09231424	Metal caps, Labocap without handle, For neck diam. 15/16 mm, Colour blue, Handle without	100 pz.
LLG09231425	Metal caps, Labocap without handle, For neck diam. 17/18 mm, Colour blue, Handle without	100 pz.
LLG09231426	Metal caps, Labocap without handle, For neck diam. 19/20 mm, Colour blue, Handle without	100 pz.
LLG09231427	Metal caps, Labocap without handle, For neck diam. 21/23 mm, Colour blue, Handle without	100 pz.
LLG09231429	Metal caps, Labocap without handle, For neck diam. 28/30 mm, Colour blue, Handle without	100 pz.
LLG09231340	Metal caps, Labocap without handle, For neck diam. 19/20 mm, Colour green, Handle without	100 pz.
LLG09231463	Metal caps, Labocap without handle, For neck diam. 12/13 mm, Colour green, Handle without	100 pz.
LLG09231464	Metal caps, Labocap without handle, For neck diam. 15/16 mm, Colour green, Handle without	100 pz.
LLG09231465	Metal caps, Labocap without handle, For neck diam. 17/18 mm, Colour green, Handle without	100 pz.



Code	Description	Packaging
LLG09231441	Metal caps, Labocap without handle, For neck diam. 9/10 mm, Colour red, Handle without	100 pz.
LLG09231442	Metal caps, Labocap without handle, For neck diam. 12/13 mm, Colour red, Handle without	100 pz.
LLG09231443	Metal caps, Labocap without handle, For neck diam. 14/15 mm, Colour red, Handle without	100 pz.
LLG09231444	Metal caps, Labocap without handle, For neck diam. 15/16 mm, Colour red, Handle without	100 pz.
LLG09231445	Metal caps, Labocap without handle, For neck diam. 17/18 mm, Colour red, Handle without	100 pz.
LLG09231446	Metal caps, Labocap without handle, For neck diam. 19/20 mm, Colour red, Handle without	100 pz.
LLG09231447	Metal caps, Labocap without handle, For neck diam. 21/23 mm, Colour red, Handle without	100 pz.
LLG09231369	Metal caps, Labocap without handle, For neck diam. 09/10 mm, Colour silver, Handle without	100 pz.
LLG09231372	Metal caps, Labocap without handle, For neck diam. 12/13 mm, Colour silver, Handle without	100 pz.
LLG09231374	Metal caps, Labocap without handle, For neck diam. 14/15 mm, Colour silver, Handle without	100 pz.
LLG09231375	Metal caps, Labocap without handle, For neck diam. 15/16 mm, Colour silver, Handle without	100 pz.
LLG09231377	Metal caps, Labocap without handle, For neck diam. 17/18 mm, Colour silver, Handle without	100 pz.
LLG09231379	Metal caps, Labocap without handle, For neck diam. 19/20 mm, Colour silver, Handle without	100 pz.
LLG09231381	Metal caps, Labocap without handle, For neck diam. 21/23 mm, Colour silver, Handle without	100 pz.
LLG09231388	Metal caps, Labocap without handle, For neck diam. 28/30 mm, Colour silver, Handle without	100 pz.
LLG09231335	Metal caps, Labocap without handle, For neck diam. 12/13 mm, Colour violet, Handle without	100 pz.
LLG09231336	Metal caps, Labocap without handle, For neck diam. 15/16 mm, Colour violet, Handle without	100 pz.
LLG09231338	Metal caps, Labocap without handle, For neck diam. 17/18 mm, Colour violet, Handle without	100 pz.
LLG09231342	Metal caps, Labocap without handle, For neck diam. 19/20 mm, Colour violet, Handle without	100 pz.
LLG09231339	Metal caps, Labocap without handle, For neck diam. 19/20 mm, Colour yellow, Handle without	100 pz.



Code	Description	Packaging
LLG09231460	Metal caps, Labocap without handle, For neck diam. 12/13 mm, Colour yellow, Handle without	100 pz.
LLG09231461	Metal caps, Labocap without handle, For neck diam. 15/16 mm, Colour yellow, Handle without	100 pz.
LLG09231462	Metal caps, Labocap without handle, For neck diam. 17/18 mm, Colour yellow, Handle without	100 pz.
LLG09231435	LABOCAP-caps 12/13mm, amberw/o hand/100	100 pz.
LLG09231438	LABOCAP caps 19/20mm, amberw/o hand/100	100 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

