

**Test tube shakers Multi Reax****Heidolph****Test tube shakers Multi Reax**

Versatile Test tube shaker with 3 mm vibration orbit, accomplishes excellent results especially in large vessels. The Multi Reax allows for rigid to smooth dispersion of viscous media and solids.

- Short term and continuous run
- Stepless speed adjustment
- Digital display
- Timer function allows for unattended operation and can be set from 1 to 999 minutes. With acoustic alarm when the operation stops.

Scope of delivery: Test tube shaker incl. 2 attachments.

- Attachment for 26 vessels with diameters from 10 mm to 16 mm
- Attachment for 12 vessels with diameters from 16 mm to 32 mm

Specifications

Type of movement: orbital

Orbit: 3 mm

Speed range:

Power consumption: 150 ... 2000 rpm

50 W

Permissible ambient temperature: 5 ... 31 °C

Permissible relative moisture: 80 %

Dimensions (W x D x H): 270 x 410 x 172 mm

Weight: 9.8 kg

Power supply: 230 ... 240 V / 50/60 Hz

IP code: IP 30



| Code | Description | Packaging |
|-------------|---|-----------|
| LLG09739259 | Accessories for Test tube shakers Multi Reax, Carousel for 12 test tubes (16 to 32 mm dia.) | 1 pz. |
| LLG09839810 | Shaker Multi Reax with rubber mat, useful area 290x258 mm | 1 pz. |
| LLG09739258 | Accessories for Test tube shakers Multi Reax, Carousel for 26 test tubes (10 to 16 mm dia.) | 1 pz. |



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

