




Hydrometers, density

Hydrometers, density

Shot poised. Without thermometer. Reference temperature: +20°C.



Code	Description	Packaging
LLG09004150	Density hydrometers without thermomet/1	
LLG09004151	Density hydrometers without thermomet/1	
LLG09004152	Density hydrometers without thermomet/1	
LLG09004153	Density hydrometers,without thermometer	1 pz.
LLG09004157	Density hydrometers,without thermometer	1 pz.
LLG09004158	Density hydrometers,without thermometer	1 pz.
LLG09004159	Density hydrometers,without thermometer	1 pz.
LLG09004161	Density hydrometers,without thermometer	1 pz.
LLG09004162	Density hydrometers,without thermometer	1 pz.
LLG09004163	Density hydrometers,without thermometer	1 pz.
LLG09004164	Density hydrometers without thermomet/1	
LLG09004165	Density hydrometers without thermomet/1	
LLG09004166	Density hydrometers without thermomet/1	
LLG09004167	Density hydrometers without thermomet/1	
LLG09004168	Density hydrometers without thermomet/1	
LLG09004170	Density hydrometers without thermomet/1	
LLG09004171	Density hydrometers without thermomet/1	
LLG09004172	Density hydrometers without thermomet/1	
LLG09004154	Density hydrometers,without thermometer	1 pz.
LLG09004155	Density hydrometers,without thermometer	1 pz.



Code	Description	Packaging
LLG09004156	Density hydrometers,without thermometer	1 pz.
LLG09004160	Density hydrometers,without thermometer	1 pz.
LLG09004169	Density hydrometers without thermomet/1	
LLG09004173	Densitiy hydrometers without thermome/1	



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

