

**Precision balances Practum**

Sartorius Lab Instruments

Precision balances Practum®

Unmatched reliability with the world-class weighing instrument. Made in Germany.

- Overload protection for years of reliability
- Supervisor lock for protection against unintentional changes
- Coloured Touchscreen-Display
- Simple data transfer to printer or PC ("Plug & Work" Technology). Direct data transfer to Microsoft® Windows programs without further software
- Integrated underfloor weighing
- Calibration: with external weights. Please order calibration weight separately
- Practum® 213/313/513: For cleaning of the windshield, the side windows, the sliding cover and all stainless steel parts are easily removable
- Extra large adjustable feet for easy levelling

Specifications


Overall dimensions (W x D x H): 360 x 216 x 95mm

PRACTUM213/313/513: 360 x 216 x 320mm

Height weighing space PRACTUM213/313/513: 209mm

Permissible ambient temperature: +10 to +30°C

Power supply: 100 - 240V, 50/60Hz



Code	Description	Packaging
LLG06262326	Precision balances Practum®, Type PRACTUM5100-1S, 5100 g	1 pz.
LLG06262840	Precision balances Practum®, Type PRACTUM6100-1S, 6100 g	1 pz.
LLG06262324	Precision balances Practum®, Type PRACTUM1102-1S, 1100 g	1 pz.
LLG06262322	Precision balances Practum®, Type PRACTUM213-1S*, 210 g	1 pz.
LLG06262323	Precision balances Practum®, Type PRACTUM2102-1S, 2100 g	1 pz.
LLG06262838	Precision balances Practum®, Type PRACTUM2101-1S, 2100 g	1 pz.
LLG06262321	Precision balances Practum®, Type PRACTUM313-1S*, 310 g	1 pz.
LLG06262837	Precision balances Practum®, Type PRACTUM3102-1S, 3100 g	1 pz.
LLG06262836	Precision balances Practum®, Type PRACTUM412-1S, 410 g	1 pz.
LLG07671426	Precision balances Practum®, Type PRACTUM5101-1S, 5100 g	1 pz.
LLG06262325	Precision balances Practum®, Type PRACTUM612-1S, 610 g	1 pz.
LLG06262839	Precision balances Practum®, Type PRACTUM6101-1S, 6100 g	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

