

**Mini mill PULVERISETTE 23****Fritsch****Mini mill PULVERISETTE 23**

Ideal for fine comminution of smallest quantities for wet and dry grinding as well as for cryogenic grinding. Its special spherical grinding bowl ensures much better performance in grinding, mixing and homogenising compared with similar models for example in chemical analysis, chromatography, mass spectrometry or x-ray structural analysis, preparation of samples for genetic research, homogenising of active ingredients, preparation of extracts in the forensic analysis, grinding of pigments and precious materials.

Max. feed size: 6mm

Max. sample quantity: up to 5ml

Final fineness: 5 μ m

Supply requirements: 100-240V/1~, 50-60Hz, 20 W

Dimensions (W x D x H): 200 x 300 x 300mm

Net weight 7kg



Code	Description	Packaging
LLG09739345	Pulverisette 23,w/o.accessories 240V/100W	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

