

**Dispensers, VITLAB Simplex2 Fix****VITLAB****Dispensers, VITLAB® Simplex2 Fix**

Fixed volume and calibration mechanism.

- with GL 45 thread
- screwable discharge valve

**Scope of supply:** 5 threaded adapters\* made of PP, 1 telescopic filling tube, instruction manual, quality certificate, 1 assembly tool.



Code	Description	Packaging
LLG06283973	Dispensers, VITLAB® Simplex2 Fix, Volume 1.0 ml, Adapter GL 28, 45, S 40	1 pz
LLG06283974	Dispensers, VITLAB® Simplex2 Fix, Volume 5.0 ml, Adapter GL 28, 45, S 40	1 pz
LLG06283975	Dispensers, VITLAB® Simplex2 Fix, Volume 10.0 ml, Adapter GL 28, 45, S 40	1 pz



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

