

**Portable Balances, Scout SKX****OHAUS****Portable Balances, Scout® SKX**

- Superior overload protection and practical stackable design for convenient storage
- Bright backlit LCD and 4 button operation make the new scout easy to see and use
- Multiple connectivity options are available for data logging and communication
- Fast weighing speed and high resolution delivers repeatable and reliable results
- Application: Weighing, density determination, mole weighing
- Display: Backlit LCD: 6-digit 7-segment with white LED backlight
- Operation: AC adapter (included), batteries (not included)
- Communication: Optional RS232, USB host, USB device, Bluetooth or Ethernet
- Construction: ABS housing, stainless steel pan, replaceable in-use cover
- Design Features: Fast stabilization time, menu lock switch, stackable design, security slot, integral weigh below hook for below balance weighing applications, transportation lock, removable stainless steel weighing platform, stability indicator, overload and underload indicators, low battery indicator, auto shut-off



Code	Description	Packaging
LLG06281345	Portable Balances, Scout® SKX, Type SKX222, 220 g, Readability 0.01 g, Plate dim. diam.120 mm	1 pz
LLG06281344	Portable Balances, Scout® SKX, Type SKX123, 120 g, Readability 0.001 g, Plate dim. diam. 93 mm	1 pz.
LLG06281349	Portable Balances, Scout® SKX, Type SKX2202, 2200 g, Readability 0.01 g, Plate dim. 170 x 140 mm	1 pz
LLG06281352	Portable Balances, Scout® SKX, Type SKX2201, 2200 g, Readability 0.1 g, Plate dim. 170 x 140 mm	1 pz
LLG06281346	Portable Balances, Scout® SKX, Type SKX422, 420 g, Readability 0.01 g, Plate dim. diam. 120 mm	1 pezzo
LLG06281350	Portable Balances, Scout® SKX, Type SKX421, 420 g, Readability 0.1 g, Plate dim. diam. 120 mm	1 pz.
LLG06281347	Portable Balances, Scout® SKX, Type SKX622, 620 g, Readability 0.01 g, Plate dim. diam. 120 mm	1 pz.
LLG06281351	Portable Balances, Scout® SKX, Type SKX621, 620 g, Readability 0.1 g, Plate dim. 170 x 140 mm	1 pz
LLG06281353	Portable Balances, Scout® SKX, Type SKX6201, 6200 g, Readability 0.1 g, Plate dim. 170 x 140 mm	1 pz.
LLG06281354	Portable Balances, Scout® SKX, Type SKX8200, 8200 g, Readability 1 g, Plate dim. 170 x 140 mm	1 pz.
LLG06281348	Portable Balances, Scout® SKX, Type SKX1202, 1200 g, Readability 0.01 g, Plate dim. 170 x 140 mm	1 pezzo



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

