

**Analytical balances Practum**

Sartorius Lab Instruments

**Analytical balances Practum®**

Unmatched reliability with the world-class weighing instrument. Made in Germany.

- Overload protection for years of reliability
- Supervisor lock for protection against unintentional changes
- Coloured Touchscreen Display
- Simple data transfer to printer or PC ("Plug & Work" Technology). Direct data transfer to Microsoft® Windows programs without further software
- Integrated underfloor weighing
- Calibration: with external weights. Please order calibration weight separately
- Easy cleaning: For cleaning of the windshield, the side windows, the sliding cover and all stainless steel parts are easily removable
- Extra large adjustable feet for easy levelling


**Specifications**

**Overall dimensions (W x D x H): 216 x 360 x 320mm**

**Height weighing space: 209mm**

**Permissible ambient temperature: +10 to +30°C**

**Power supply: 100 - 240V, 50/60Hz**



Code	Description	Packaging
LLG06260896	Analytical balances Practum®, Type PRACTUM224-1S, Weighing range 220 g	1 pz.
LLG06262320	Analytical balances Practum®, Type PRACTUM124-1S, Weighing range 120 g	1 pz.
LLG06262847	Analytical balances Practum®, Type PRACTUM64-1S, Weighing range 60 g	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

