



## Accessories for collecting pans, stainless steel

### Accessories for collecting pans, stainless steel

Stainless steel. In accordance to ISO 3310-1:2016 and ASTM E11:2009.



Code	Description	Packaging
LLG09739237	Accessories for collecting pans, stainless steel, 100mm dia., with outlet	1 pz.
LLG09739236	Accessories for collecting pans, stainless steel, 200mm dia., with outlet	1 pz.
LLG09739235	Accessories for collecting pans, stainless steel, 200mm dia., without outlet	1 pz.
LLG09739240	Accessories for collecting pans, stainless steel, 203mm dia., with outlet	1 pz.
LLG09739239	Accessories for collecting pans, stainless steel, 203mm dia., without outlet	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

