



Portable Balances, Scout STX

OHAUS

Portable Balances, Scout® STX

- Navigate quickly and efficiently using the informative color touchscreen display
- Fast weighing speed and high resolution deliver repeatable and reliable results
- Superior overload protection and stackable storage drive overall durability
- Multiple weighing modes and languages for simplified weighing
- Application: Weighing, percent weighing, parts counting, check weighing, dynamic/animal weighing, totalization, density determination, display hold, mole weighing
- Display: Full-color touchscreen WQVGA graphic LCD (4.3" diagonal), QWERTY keyboard and numeric keypad to quickly input GLP and GMP data and other application data
- Operation: AC adapter (Included), batteries (not included)
- Communication: Optional RS232, USB host, USB device, Bluetooth or Ethernet
- Construction: ABS housing, stainless steel pan, replaceable in-use cover
- Design Features: Fast stabilization time, up-front level indicator, up to 9 operating languages, save-to-USB function, menu lock switch, stackable design, security slot, integral weigh below hook for below balance weighing applications, transportation lock, removable stainless steel weighing platform, stability indicator, overload and underload indicators, low battery indicator, auto shut-off



Code	Description	Packaging
LLG06281360	Portable Balances, Scout® STX, Type STX1202, 1200 g, Readability 0.01 g, Plate dim. 170 x 140 mm	1 pz.
LLG06281356	Portable Balances, Scout® STX, Type STX223, 220 g, Readability 0.001 g, Plate dim. diam. 93 mm	1 pz.
LLG06281357	Portable Balances, Scout® STX, Type STX222, 220 g, Readability 0.01 g, Plate dim. diam. 120 mm	1 pz.
LLG06281361	Portable Balances, Scout® STX, Type STX2202, 2200 g, Readability 0.01 g, Plate dim. 170 x 140 mm	1 pz.
LLG06281364	Portable Balances, Scout® STX, Type STX2201, 2200 g, Readability 0.1 g, Plate dim. 170 x 140 mm	1 pz.
LLG06281358	Portable Balances, Scout® STX, Type STX422, 420 g, Readability 0.01 g, Plate dim. diam. 120 mm	1 pz.
LLG06281362	Portable Balances, Scout® STX, Type STX421, 420 g, Readability 0.1 g, Plate dim. diam. 120 mm	1 pz.
LLG06281359	Portable Balances, Scout® STX, Type STX622, 620 g, Readability 0.01 g, Plate dim. diam. 120 mm	1 pz.
LLG06281363	Portable Balances, Scout® STX, Type STX621, 620 g, Readability 0.1 g, Plate dim. 170 x 140 mm	1 pz.
LLG06281366	Portable Balances, Scout® STX, Type STX8200, 8200 g, Readability 1 g, Plate dim. 170 x 140 mm	1 pz.
LLG06281355	Portable Balances, Scout® STX, Type STX123, 120 g, Readability 0.001 g, Plate dim. diam. 93 mm	1 pz.
LLG06281365	Portable Balances, Scout® STX, Type STX6201, 6200 g, Readability 0.1 g, Plate dim. 170 x 140 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

