

**Pocket Balances Type JE****OHAUS****Pocket Balances Type JE**

Simple, compact and handy. Ideal for precise weighing in the field.

- Automatic calibration, using external weight
- Large, high-contrast, liquid crystal display with back light
- Robust plastic housing
- Auto-off function prolongs battery life
- Includes protective cover and batteries
- Dimensions (W x D x H): 136 x 83 x 20mm



Code	Description	Packaging
LLG09902413	Pocket Balances Type JE, Type JE120, 120 g, Readability 0.1 g, Plate dim. 83 x 76 mm	1 pz.
LLG09902414	Pocket Balances Type JE, Type JE250, 250 g, Readability 0.1 g, Plate dim. 83 x 76 mm	1 pz.
LLG09902415	Pocket Balances Type JE, Type JE500, 500 g, Readability 0.1 g, Plate dim. 83 x 76 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

