



Deep-well plate system, Riplate

Deep-well plate system, Riplate®

Riplate® round wells - Standard in Multifunction

Raised wells reduce the risk of cross-contamination and increase the security of hot sealing or sealing with adhesive foils. Riplate® 1ml is ideal for use in a thermocycler or for defrosting of samples in a hot-water bath. The U-shaped well bottom reduces the death-volume.

Riplate® square wells - For a Maximum in Capacity

Riplate® SW have been developed for optimal use of capacity when it comes to sample storage. Due to the SBS-format the outside dimensions of Riplate® are given. Form follows function: the construction of cavities is crucial for a maximum absorption of liquids.

Riplate® medio - Saves 50% of Storage Capacity

Due to the compact construction of Riplate® medio less than 50% stacking size is needed. Higher well rims allow an airtight closure or sealing of Riplate® medio. The alphanumeric marking allows registration and identification of single assays.



Riplate® 384 - Optimized for Smaller Volumes

Available in PS and PP for automated screening and storage systems. Alphanumeric marking for identification of samples. Round cavities for reduced capillary forces. With flat well-bottom - ideal for microscopic and optical applications.

Field of application: screening, sample storage et al. Ritter Riplate® systems are in line with the SBS-standard and are suitable for automation systems with corresponding closing- and sealing-mechanism.



Code	Description	Packaging
LLG09284228	Deep-well plate system, Riplate®, Type Cap mat for Square well 48, 5ml	100 pz.
LLG09284209	Deep-well plate system, Riplate®, Type Cap mat for Square well 96, 2ml	100 pz.
LLG09284296	Deep-well plate system, Riplate®, Type Deep-well 96, PP - 0.5ml medio	20 pz.
LLG09284295	Deep-well plate system, Riplate®, Type Deep-well 96, PP - 1.0ml medio	32 pz.
LLG09284251	Deep-well plate system, Riplate®, Type Round well 96, PP - 1ml, U-shaped well-bottom	20 pz.
LLG09284300	Deep-well plate system, Riplate®, Type Round well 96, PP - 2ml, U-shaped well-bottom	20 pz.
LLG09284250	Deep-well plate system, Riplate®, Type Round well 96, PS - 1ml, U-shaped well-bottom	20 pz.
LLG07635031	Deep-well plate system, Riplate®, Type Square-well 24, PP - 10ml, V-shaped well-bottom	20 pz.
LLG09284260	Deep-well plate system, Riplate®, Type Square-well 48, PP - 5ml, V-shaped well-bottom	20 pz.
LLG09284252	Deep-well plate system, Riplate®, Type Square-well 96, PP - 2ml, U-shaped well-bottom	20 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

