


**Gas sampling tubes, DURAN tubing****Lenz****Gas sampling tubes, DURAN® tubing**

Made of DURAN® tubing. Available with or without nozzle (with thread GL and septa), with two one-way stopcocks NS14 (with retaining device), connections: capillary tube (outer dia. 8mm, inner dia. 2.7mm), glass olive (outer dia. 8mm), with inscription label. Also available with PUR coating.



Code	Description	Packaging
LLG07623930	Gas sampling tubes, DURAN® tubing, Type With sampling connections, GL 14, 1000 ml, External body diam. 85 mm, Length* 260	1 pz.
LLG09110230	Gas sampling tubes, DURAN® tubing, Type Without sampling connections, 1000 ml, External body diam. 85 mm, Length* 260	1 pz.
LLG06243069	Gas sampling tubes, DURAN® tubing, Type With sampling connections, GL 14, 150 ml, External body diam. 42 mm, Length* 155	1 pz.
LLG09110215	Gas sampling tubes, DURAN® tubing, Type Without sampling connections, 150 ml, External body diam. 42 mm, Length* 155	1 pz.
LLG09110220	Gas sampling tubes, DURAN® tubing, Type Without sampling connections, 350 ml, External body diam. 54 mm, Length* 200	1 pz.
LLG09110225	Gas sampling tubes, DURAN® tubing, Type Without sampling connections, 500 ml, External body diam. 65 mm, Length* 220	1 pz.
LLG09110237	Gas sampling tubes, DURAN® tubing, Type With sampling connections, GL 14, 500 ml, External body diam. 65 mm, Length* 220	1 pz.
LLG09110236	Gas sampling tubes, DURAN® tubing, Type With sampling connections, GL 14, 350 ml, External body diam. 54 mm, Length* 200	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

