

**Beer Force-Ageing-Test Bath BFT5****Huber****Beer Force-Ageing-Test Bath BFT®5**

We offer special air cooled circulator for the Beer Force-Ageing-Test. It is equipped with a programmer for the change between 0 °C / 40 °C and 0 °C / 60 °C in a 24 hour cycle. Casings and bath parts are made of stainless steel.

Dimensions (W x D x H): 460 x 710 x 911mm



Code	Description	Packaging
LLG09857399	Beer Force-Ageing-Test Bath BFT®5, Bath opening 350 x 410 mm, Bath tank depth 270 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

