



Hoses

Hoses



Code	Description	Packaging
LLG09857046	Hoses, Type NW 12 M24 x 1.5, Temp.working range -120 ... 400 °C, Length 100 cm, Material Metal	1 pz.
LLG09857121	Hoses, Type NW 12 M16 x 1, Temp.working range -100 ... 350 °C, Length 100 cm, Material Metal	1 pz.
LLG09857125	Hoses, Type NW 20 M30 x 1.5, Temp.working range -100 ... 350 °C, Length 100 cm, Material Metal	1 pz.
LLG09857129	Hoses, Type NW 25 M38 x 1.5, Temp.working range -100 ... 350 °C, Length 100 cm, Material Metal	1 pz.
LLG09857214	Hoses, Type NW 12 M16 x 1, Temp.working range -50 ... 200 °C, Length 100 cm, Material Metal	1 pz.
LLG09857218	Hoses, Type NW 12 M24 x 1.5, Temp.working range -100 ... 350 °C, Length 100 cm, Material Metal	1 pz.
LLG09857047	Hoses, Type NW 12 M24 x 1.5, Temp.working range -120 ... 400 °C, Length 150 cm, Material Metal	1 pz.
LLG09857122	Hoses, Type NW 12 M16 x 1, Temp.working range -100 ... 350 °C, Length 150 cm, Material Metal	1 pz.
LLG09857126	Hoses, Type NW 20 M30 x 1.5, Temp.working range -100 ... 350 °C, Length 150 cm, Material Metal	1 pz.
LLG09857130	Hoses, Type NW 25 M38 x 1.5, Temp.working range -100 ... 350 °C, Length 150 cm, Material Metal	1 pz.
LLG09857215	Hoses, Type NW 12 M16 x 1, Temp.working range -50 ... 200 °C, Length 150 cm, Material Metal	1 pz.
LLG09857045	Hoses, Type NW 12 M16 x 1, Temp.working range -100 ... 350 °C, Length 200 cm, Material Metal	1 pz.
LLG09857048	Hoses, Type NW 12 M24 x 1.5, Temp.working range -120 ... 400 °C, Length 200 cm, Material Metal	1 pz.
LLG09857127	Hoses, Type NW 20 M30 x 1.5, Temp.working range -100 ... 350 °C, Length 200 cm, Material Metal	1 pz.
LLG09857131	Hoses, Type NW 25 M38 x 1.5, Temp.working range -100 ... 350 °C, Length 200 cm, Material Metal	1 pz.
LLG09857216	Hoses, Type NW 12 M16 x 1, Temp.working range -50 ... 200 °C, Length 200 cm, Material Metal	1 pz.
LLG09857219	Hoses, Type NW 12 M24 x 1.5, Temp.working range -100 ... 350 °C, Length 200 cm, Material Metal	1 pz.
LLG09857049	Hoses, Type NW 12 M24 x 1.5, Temp.working range -120 ... 400 °C, Length 300 cm, Material Metal	1 pz.
LLG09857124	Hoses, Type NW 12 M16 x 1, Temp.working range -100 ... 350 °C, Length 300 cm, Material Metal	1 pz.
LLG09857128	Hoses, Type NW 20 M30 x 1.5, Temp.working range -100 ... 350 °C, Length 300 cm, Material Metal	1 pz.



Code	Description	Packaging
LLG09857132	Hoses, Type NW 25 M38 x 1.5, Temp.working range -100 ... 350 °C, Length 300 cm, Material Metal	1 pz.
LLG09857217	Hoses, Type NW 12 M16 x 1, Temp.working range -50 ... 200 °C, Length 300 cm, Material Metal	1 pz.
LLG09857220	Hoses, Type NW 12 M24 x 1.5, Temp.working range -100 ... 350 °C, Length 300 cm, Material Metal	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

