

**LLG-Erlenmeyer flasks with standard ground joint, borosilicate glass 3.3****LLG Labware****LLG-Erlenmeyer flasks with standard ground joint, borosilicate glass 3.3**

Ideal for mixing liquids. With standard ground joint and white graduation. Acc. to DIN EN ISO 4797.

Nominal capacity 50 ml

Neck Ø

Bottom Ø 51 mm

Height 85 mm

Thread

Cap size

Ground size NS14/23

Standard DIN EN ISO 4797

Graduation

Colour clear

Material Boro 3.3

Autoclavable Yes



| Code        | Description                                                                                                                                                                                                                                                 | Packaging |
|-------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------|
| LLG04686106 | LLG-Erlenmeyer flask 50 ml pack of 2                                                                                                                                                                                                                        | 2 pz      |
| LLG04686107 | LLG-Erlenmeyer flask 50 ml, pack of 2                                                                                                                                                                                                                       | 2 pz      |
| LLG04686108 | LLG-Erlenmeyer flask 50 ml, pack of 2                                                                                                                                                                                                                       | 2 pz      |
| LLG04686109 | LLG-Erlenmeyer flask 100 ml pack of 2                                                                                                                                                                                                                       | 2 pz      |
| LLG04686110 | LLG-Erlenmeyer flasks with standard ground joint, borosilicate glass 3.3, Nominal capacity 100 ml, Neck Ø , Bottom Ø 64 mm, Height 100 mm, Thread , Cap size , Ground size NS19/26 , Standard DIN EN ISO 4797 , Graduation , Colour clear , Material Boro 3 | 2 pz      |
| LLG04686111 | LLG-Erlenmeyer flask 100 ml NS29/32 x2                                                                                                                                                                                                                      | 2 pz      |
| LLG04686112 | LLG-Erlenmeyer flasks with standard ground joint, borosilicate glass 3.3, Nominal capacity 250 ml, Neck Ø , Bottom Ø 85 mm, Height 140 mm, Thread , Cap size , Ground size NS19/26 , Standard DIN EN ISO 4797 , Graduation , Colour clear , Material Boro 3 | 2 pz      |
| LLG04686113 | LLG-Erlenmeyer flasks with standard ground joint, borosilicate glass 3.3, Nominal capacity 250 ml, Neck Ø , Bottom Ø 85 mm, Height 140 mm, Thread , Cap size , Ground size NS24/29 , Standard DIN EN ISO 4797 , Graduation , Colour clear , Material Boro 3 | 2 pz      |
| LLG04686114 | LLG-Erlenmeyer flasks with standard ground joint, borosilicate glass 3.3, Nominal capacity 250 ml, Neck Ø , Bottom Ø 85 mm, Height 140 mm, Thread , Cap size , Ground size NS29/32 , Standard DIN EN ISO 4797 , Graduation , Colour clear , Material Boro 3 | 2 pz      |
| LLG04686115 | LLG-Erlenmeyer flasks with standard ground joint, borosilicate glass 3.3, Nominal capacity 300 ml, Neck Ø , Bottom Ø 87 mm, Height 156 mm, Thread , Cap size , Ground size NS29/32 , Standard , Graduation , Colour clear , Material Boro 3.3 , Autoclavab  | 2 pz      |
| LLG04686116 | LLG-Erlenmeyer flasks with standard ground joint, borosilicate glass 3.3, Nominal capacity 500 ml, Neck Ø , Bottom Ø 105 mm, Height 175 mm, Thread , Cap size , Ground size NS29/32 , Standard DIN EN ISO 4797 , Graduation , Colour clear , Material Boro  | 2 pz      |

| Code        | Description                                                                                                                                                                                                                                                 | Packaging |
|-------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------|
| LLG04686117 | LLG-Erlenmeyer flasks with standard ground joint, borosilicate glass 3.3, Nominal capacity 1000 ml, Neck Ø , Bottom Ø 131 mm, Height 220 mm, Thread , Cap size , Ground size NS29/32 , Standard DIN EN ISO 4797 , Graduation , Colour clear , Material Boro | 1 pz      |
| LLG04686118 | LLG-Erlenmeyer flasks with standard ground joint, borosilicate glass 3.3, Nominal capacity 2000 ml, Neck Ø , Bottom Ø 166 mm, Height 270 mm, Thread , Cap size , Ground size NS29/32 , Standard DIN EN ISO 4797 , Graduation , Colour clear , Material Boro | 1 pz      |



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

