

**LLG-Erlenmeyer flasks with standard ground joint, borosilicate glass 3.3****LLG Labware****LLG-Erlenmeyer flasks with standard ground joint, borosilicate glass 3.3**

Ideal for mixing liquids. With standard ground joint and white graduation. Acc. to DIN EN ISO 4797.

Nominal capacity 50 ml
Neck Ø
Bottom Ø 51 mm
Height 85 mm
Thread
Cap size
Ground size NS14/23
Standard DIN EN ISO 4797
Graduation
Colour clear
Material Boro 3.3
Autoclavable Yes



Code	Description	Packaging
LLG04686106	LLG-Erlenmeyer flask 50 ml pack of 2	2 pz
LLG04686107	LLG-Erlenmeyer flask 50 ml, pack of 2	2 pz
LLG04686108	LLG-Erlenmeyer flask 50 ml, pack of 2	2 pz
LLG04686109	LLG-Erlenmeyer flask 100 ml pack of 2	2 pz
LLG04686110	LLG-Erlenmeyer flasks with standard ground joint, borosilicate glass 3.3, Nominal capacity 100 ml, Neck Ø , Bottom Ø 64 mm, Height 100 mm, Thread , Cap size , Ground size NS19/26 , Standard DIN EN ISO 4797 , Graduation , Colour clear , Material Boro 3	2 pz
LLG04686111	LLG-Erlenmeyer flask 100 ml NS29/32 x2	2 pz
LLG04686112	LLG-Erlenmeyer flasks with standard ground joint, borosilicate glass 3.3, Nominal capacity 250 ml, Neck Ø , Bottom Ø 85 mm, Height 140 mm, Thread , Cap size , Ground size NS19/26 , Standard DIN EN ISO 4797 , Graduation , Colour clear , Material Boro 3	2 pz
LLG04686113	LLG-Erlenmeyer flasks with standard ground joint, borosilicate glass 3.3, Nominal capacity 250 ml, Neck Ø , Bottom Ø 85 mm, Height 140 mm, Thread , Cap size , Ground size NS24/29 , Standard DIN EN ISO 4797 , Graduation , Colour clear , Material Boro 3	2 pz
LLG04686114	LLG-Erlenmeyer flasks with standard ground joint, borosilicate glass 3.3, Nominal capacity 250 ml, Neck Ø , Bottom Ø 85 mm, Height 140 mm, Thread , Cap size , Ground size NS29/32 , Standard DIN EN ISO 4797 , Graduation , Colour clear , Material Boro 3	2 pz
LLG04686115	LLG-Erlenmeyer flasks with standard ground joint, borosilicate glass 3.3, Nominal capacity 300 ml, Neck Ø , Bottom Ø 87 mm, Height 156 mm, Thread , Cap size , Ground size NS29/32 , Standard , Graduation , Colour clear , Material Boro 3.3 , Autoclavab	2 pz
LLG04686116	LLG-Erlenmeyer flasks with standard ground joint, borosilicate glass 3.3, Nominal capacity 500 ml, Neck Ø , Bottom Ø 105 mm, Height 175 mm, Thread , Cap size , Ground size NS29/32 , Standard DIN EN ISO 4797 , Graduation , Colour clear , Material Boro	2 pz



Code	Description	Packaging
LLG04686117	LLG-Erlenmeyer flasks with standard ground joint, borosilicate glass 3.3, Nominal capacity 1000 ml, Neck Ø , Bottom Ø 131 mm, Height 220 mm, Thread , Cap size , Ground size NS29/32 , Standard DIN EN ISO 4797 , Graduation , Colour clear , Material Boro	1 pz
LLG04686118	LLG-Erlenmeyer flasks with standard ground joint, borosilicate glass 3.3, Nominal capacity 2000 ml, Neck Ø , Bottom Ø 166 mm, Height 270 mm, Thread , Cap size , Ground size NS29/32 , Standard DIN EN ISO 4797 , Graduation , Colour clear , Material Boro	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

