

**LLG-Erlenmeyer flasks, borosilicate glass 3.3, with screw cap****LLG Labware****LLG-Erlenmeyer flasks, borosilicate glass 3.3, with screw cap**

Erlenmeyer flasks made of borosilicate glass 3.3. Screw cap with PTFE liner.

- White graduation and labelling field
- GL thread

Nominal capacity 250 ml

Neck Ø

Bottom Ø 85 mm

Height 140 mm

Thread GL32

Cap size

Ground size

Standard

Graduation

Colour clear

Material Boro 3.3

Autoclavable Yes



Code	Description	Packaging
LLG04686198	LLG-Erlenmeyer flasks, borosilicate glass 3.3, with screw cap, Nominal capacity 100 ml, Neck Ø , Bottom Ø 64 mm, Height 100 mm, Thread GL32 , Cap size , Ground size , Standard , Graduation , Colour clear , Material Boro 3.3 , Autoclavable Yes	2 pz
LLG04686199	LLG-Erlenmeyer flasks, borosilicate glass 3.3, with screw cap, Nominal capacity 250 ml, Neck Ø , Bottom Ø 85 mm, Height 140 mm, Thread GL32 , Cap size , Ground size , Standard , Graduation , Colour clear , Material Boro 3.3 , Autoclavable Yes	2 pz
LLG04686200	LLG-Erlenmeyer flask 500 ml, boro 3.3, w	2 pz
LLG04686201	LLG-Erlenmeyer flask 1000 ml, boro 3.3,	1 pz



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

