

**LLG-Measuring cylinders, borosilicate glass 3.3, tall form, class A****LLG Labware****LLG-Measuring cylinders, borosilicate glass 3.3, tall form, class A**

- Tall form
- Conformity certified with batch ID.
- DIN EN ISO 4788
- With hexagonal base and spout
- Blue graduation



Code	Description	Packaging
LLG04686202	LLG-Measuring cylinder 5 ml, pack 2	2 pz
LLG04686203	LLG-Measuring cylinder 10 ml, pack 2	2 pz
LLG04686204	LLG-Measuring cylinder 25 ml, pack 2	2 pz
LLG04686205	LLG-Measuring cylinder 50 ml, pack 2	2 pz
LLG04686206	LLG-Measuring cylinder 100 ml, pack 2	2 pz
LLG04686207	LLG-Measuring cylinder 250 ml, pack 2	2 pz
LLG04686208	LLG-Measuring cylinder 500 ml, pack 2	2 pz
LLG04686209	LLG-Measuring cylinder 1000 ml,	1pz
LLG04686210	LLG-Measuring cylinder 2000 ml,	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

