

**LLG-Measuring cylinders, borosilicate glass 3.3, low form, class B****LLG Labware****LLG-Measuring cylinders, borosilicate glass 3.3, low form, class B**

- Low form
- Calibrated to contain (TC, IN)
- DIN EN ISO 4788
- With hexagonal base and spout
- Amber graduation



Code	Description	Packaging
LLG04686211	LLG-Measuring cylinders 10ml, pack 2	2 pz
LLG04686212	LLG-Measuring cylinders 25ml, pack 2	2 pz
LLG04686213	LLG-Measuring cylinders 25ml, pack 2	2 pz
LLG04686214	LLG-Measuring cylinders 100ml, pack 2	2 pz
LLG04686215	LLG-Measuring cylinders 250ml, pack 2	2 pz
LLG04686216	LLG-Measuring cylinders 500ml, pack 2	2 pz
LLG04686217	LLG-Measuring cylinder 1000ml	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

