



**LLG-Mixing cylinders, borosilicate glass 3.3, tall form, class A**

**LLG Labware**

### LLG-Mixing cylinders, borosilicate glass 3.3, tall form, class A

- Tall form
- Conformity certified with batch ID
- DIN EN ISO 4788
- With hexagonal base and PE stopper
- Blue graduation



Code	Description	Packaging
LLG04686218	LLG-Mixing cylinders 10 ml, pack of 2	2 pz
LLG04686219	LLG-Mixing cylinders 25 ml, pack of 2	2 pz
LLG04686220	LLG-Mixing cylinders 50 ml, pack of 2	2 pz
LLG04686221	LLG-Mixing cylinders 100 ml, pack of 2	2 pz
LLG04686222	LLG-Mixing cylinders 250 ml, pack of 2	2 pz
LLG04686223	LLG-Mixing cylinders 500 ml, pack of 2	2 pz
LLG04686224	LLG-Mixing cylinders 1000 ml, pack of 2	1 pz



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

