

**LLG-Volumetric flasks, borosilicate glass 3.3, class A****LLG Labware****LLG-Volumetric flasks, borosilicate glass 3.3, class A**

Conformity certified with batch ID.

- DIN EN ISO 1042
- With standard ground joint and PE stopper
- Blue ring mark and inscription

Nominal capacity 25 ml

Tolerance ± 0.040 ml

Ground size NS10/19



| Code | Description | Packaging |
|-------------|---|-----------|
| LLG04686232 | LLG-Volumetric flask 5 ml, pack of 2 | 2 pz |
| LLG04686233 | LLG-Volumetric flask 10 ml, pack of 2 | 2 pz |
| LLG04686234 | LLG-Volumetric flask 20 ml, pack of 2 | 2 pz |
| LLG04686235 | LLG-Volumetric flask 25 ml, pack of 2 | 2 pz |
| LLG04686236 | LLG-Volumetric flask 25 ml, pack of 2 | 2 pz |
| LLG04686237 | LLG-Volumetric flask 50 ml, pack of 2 | 2 pz |
| LLG04686238 | LLG-Volumetric flask 50 ml, pack of 2 | 2 pz |
| LLG04686239 | LLG-Volumetric flask 100 ml, pack of 2 | 2 pz |
| LLG04686240 | LLG-Volumetric flask 100 ml, pack of 2 | 2 pz |
| LLG04686241 | LLG-Volumetric flask 200 ml, pack of 2 | 2 pz |
| LLG04686242 | LLG-Volumetric flask 250 ml, pack of 2 | 2 pz |
| LLG04686243 | LLG-Volumetric flask 500 ml, pack of 2 | 2 pz |
| LLG04686244 | LLG-Volumetric flask 1000 ml, pack of 2 | 1 pz |
| LLG04686245 | LLG-Volumetric flask 2000 ml, pack of 2 | 1 pz |





CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

