

**LLG-Volumetric flasks, borosilicate glass 3.3, class A****LLG Labware****LLG-Volumetric flasks, borosilicate glass 3.3, class A**

Conformity certified with batch ID.

- DIN EN ISO 1042
- With standard ground joint and PE stopper
- Blue ring mark and inscription

Nominal capacity 25 ml

Tolerance  $\pm 0.040$  ml

Ground size NS10/19



Code	Description	Packaging
LLG04686232	LLG-Volumetric flask 5 ml, pack of 2	2 pz
LLG04686233	LLG-Volumetric flask 10 ml, pack of 2	2 pz
LLG04686234	LLG-Volumetric flask 20 ml, pack of 2	2 pz
LLG04686235	LLG-Volumetric flask 25 ml, pack of 2	2 pz
LLG04686236	LLG-Volumetric flask 25 ml, pack of 2	2 pz
LLG04686237	LLG-Volumetric flask 50 ml, pack of 2	2 pz
LLG04686238	LLG-Volumetric flask 50 ml, pack of 2	2 pz
LLG04686239	LLG-Volumetric flask 100 ml, pack of 2	2 pz
LLG04686240	LLG-Volumetric flask 100 ml, pack of 2	2 pz
LLG04686241	LLG-Volumetric flask 200 ml, pack of 2	2 pz
LLG04686242	LLG-Volumetric flask 250 ml, pack of 2	2 pz
LLG04686243	LLG-Volumetric flask 500 ml, pack of 2	2 pz
LLG04686244	LLG-Volumetric flask 1000 ml, pack of 2	1 pz
LLG04686245	LLG-Volumetric flask 2000 ml, pack of 2	1 pz





**CARLO ERBA Reagents S.A.S.**  
Chaussée du Vexin – 27106 Val  
de Reuil cedex  
N° TVA: FR 63391048824  
Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

