



LLG-Volumetric flasks, borosilicate glass 3.3, class A

LLG Labware

LLG-Volumetric flasks, borosilicate glass 3.3, class A

Conformity certified with batch ID.

- DIN EN ISO 1042
- With standard ground joint and PE stopper
- Blue ring mark and inscription

Nominal capacity 25 ml

Tolerance ± 0.040 ml

Ground size NS10/19



Code	Description	Packaging
LLG04686232	LLG-Volumetric flask 5 ml, pack of 2	2 pz
LLG04686233	LLG-Volumetric flask 10 ml, pack of 2	2 pz
LLG04686234	LLG-Volumetric flask 20 ml, pack of 2	2 pz
LLG04686235	LLG-Volumetric flask 25 ml, pack of 2	2 pz
LLG04686236	LLG-Volumetric flask 25 ml, pack of 2	2 pz
LLG04686237	LLG-Volumetric flask 50 ml, pack of 2	2 pz
LLG04686238	LLG-Volumetric flask 50 ml, pack of 2	2 pz
LLG04686239	LLG-Volumetric flask 100 ml, pack of 2	2 pz
LLG04686240	LLG-Volumetric flask 100 ml, pack of 2	2 pz
LLG04686241	LLG-Volumetric flask 200 ml, pack of 2	2 pz
LLG04686242	LLG-Volumetric flask 250 ml, pack of 2	2 pz
LLG04686243	LLG-Volumetric flask 500 ml, pack of 2	2 pz
LLG04686244	LLG-Volumetric flask 1000 ml, pack of 2	1 pz
LLG04686245	LLG-Volumetric flask 2000 ml, pack of 2	1 pz





CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

