



LLG-Volumetric flasks, borosilicate glass 3.3, class A, amber glass

LLG Labware

LLG-Volumetric flasks, borosilicate glass 3.3, class A, amber glass

- Conformity certified with batch ID
- DIN EN ISO 1042
- With standard ground joint and PE topper
- White ring mark and inscription



| Code | Description | Packaging |
|-------------|---|-----------|
| LLG04686246 | LLG-Volumetric flasks 5 ml,, pack of 2 | 2 pz |
| LLG04686247 | LLG-Volumetric flasks 10 ml,, pack of 2 | 2 pz |
| LLG04686248 | LLG-Volumetric flasks 20 ml,, pack of 2 | 2 pz |
| LLG04686249 | LLG-Volumetric flasks 25 ml,, pack of 2 | 2 pz |
| LLG04686250 | LLG-Volumetric flasks 50 ml,, pack of 2 | 2 pz |
| LLG04686251 | LLG-Volumetric flasks 100 ml pack of 2 | 2 pz |
| LLG04686252 | LLG-Volumetric flasks 200 ml pack of 2 | 2 pz |
| LLG04686253 | LLG-Volumetric flask 250 ml, pack of 2 | 2 pz |
| LLG04686254 | LLG-Volumetric flask 500 ml, pack of 2 | 2 pz |
| LLG04686255 | LLG-Volumetric flask 1000 ml, pack of 2 | 1 pz |



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

