

**LLG-Volumetric trapezoidal flasks, borosilicate glass 3.3, class A****LLG Labware****LLG-Volumetric trapezoidal flasks, borosilicate glass 3.3, class A**

- Trapezoidal shape
- Conformity certified with batch ID
- DIN EN ISO 1042
- With standard ground joint and PE stopper
- Blue ring mark and inscription



Code	Description	Packaging
LLG04686225	LLG-Volumetric flask 1 ml, pack of 2	2 pz
LLG04686226	LLG-Volumetric flask 2 ml, pack of 2	2 pz
LLG04686227	LLG-Volumetric flask 3 ml, pack of 2	2 pz
LLG04686228	LLG-Volumetric flask 4 ml, pack of 2	2 pz
LLG04686229	LLG-Volumetric flask 10 ml, pack of 2	2 pz
LLG04686230	LLG-Volumetric flask 20 ml, pack of 2	2 pz
LLG04686231	LLG-Volumetric flask 25 ml, pack of 2	2 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

