

**LLG-Separating funnel acc. to Squibb, borosilicate glass 3.3****LLG Labware****LLG-Separating funnel acc. to Squibb, borosilicate glass 3.3**

Separating funnel acc. to Squibb with PTFE stopcock.

- DIN ISO 4800
- With standard ground joint and stopper
- White graduation



Code	Description	Packaging
LLG04686263	LLG-Separating funnel 50 ml Squibb	1 pz
LLG04686264	LLG-Separating funnels 100 ml Squibb	1 pz
LLG04686265	LLG-Separating funnels 250 ml Squibb	1 pz
LLG04686266	LLG-Separating funnels 500 ml Squibb	1 pz
LLG04686267	LLG-Separating funnels 1000 ml Squibb	1 pz
LLG04686268	LLG-Separating funnels 2000 ml Squibb	1 pz



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

