

**LLG-Hexagonal hollow stoppers, borosilicate glass 3.3, pointed****LLG Labware****LLG-Hexagonal hollow stoppers, borosilicate glass 3.3, pointed**

For closing glassware with standard ground joints. The drip tip ensures that adhering liquid flows back into the glassware. Hexagonal shape prevents the stopper from rolling away. Acc. to DIN 12252.

Ground size NS14/23



Code	Description	Packaging
LLG04686100	LLG-Hexagonal stopper 10/19, pack of 5	5 pz
LLG04686101	LLG-Hexagonal stopper 10/19, pack of 5	5 pz
LLG04686102	LLG-Hexagonal stopper NS 14/23, boro 3.3	5 pz
LLG04686103	LLG-Hexagonal stopper 19/26, pack of 5	5 pz
LLG04686104	LLG-Hexagonal hollow stoppers, borosilicate glass 3.3, pointed, Ground size NS24/29	5 pz
LLG04686105	LLG-Hexagonal hollow stoppers, borosilicate glass 3.3, pointed, Ground size NS29/32	5 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

