

**LLG-Condenser acc. to Dimroth, borosilicate glass 3.3****LLG Labware****LLG-Condenser acc. to Dimroth, borosilicate glass 3.3, glass olive or PP side**

- For condensing vapours.
- Water connection with glass olives or PP side
- With standard ground cone and socket
- Acc. to DIN 12591



Code	Description	Packaging
LLG04686190	LLG-Dimroth condenser 250 mm, glass olive	1 pz
LLG04686191	LLG-Dimroth condenser 400 mm, glass olive	1 pz
LLG04686192	LLG-Dimroth condenser 250mm PP side arm	1 pz
LLG04686193	LLG-Dimroth condenser 400mm PP side arm	1 pz



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

