

**LLG-Condenser acc. to Liebig, borosilicate glass 3.3, glass olive****LLG Labware****LLG-Condenser acc. to Liebig, borosilicate glass 3.3, glass olive**

For condensing vapours. Water connection with glass olives. With standard ground cone and socket. Acc. to DIN 12576.



Code	Description	Packaging
LLG04686188	LLG-Liebig condenser 250mm NS29/32	1 pz
LLG04686189	LLG-Liebig condenser 400mm 29/32 glass	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

