



LLG-Crimp seals ND20, center tear off caps, aluminium, unassembled

LLG Labware

LLG-Crimp seals ND20, center tear off caps, aluminium, unassembled



Code	Description	Packaging
LLG06270720	LLG-Crimp seals ND20, center tear off caps, aluminium, unassembled, Colour Silver, With hole Yes, Int. bore Ø 8.0 mm, Material septa Butyl, Colour septa Grey, Thickness 3.00 mm, Shore hardness 50 A	100 pz



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

