

**Magnetic stirrers C-MAG MS 4 / C-MAG MS 7 / C-MAG MS 10****Magnetic stirrers C-MAG MS 4 / C-MAG MS 7 / C-MAG MS 10**

Without heating. With glass ceramic top plate which offers excellent chemical resistance.

-powerful motor for stirring quantities of up to 5l, 10l, or 15l (H₂O)

-raised control panel for protection against leaking liquids

Specifications**Magnetic stirrer**

Stirring volume max. (H₂O)

C-Mag MS 4 / MS 7 / MS 10:

5l / 10l / 15l

Motor rating input / output: 15 / 1.5W

Speed range: 100 to 1500rpm

Speed display:

Max. magnetic bar (L x dia.)

C-MAG MS 4:

C-MAG MS 7 / MS 10: scale

30 x 8mm

80 x 10mm

Set-up plate

Material:

glass ceramic

Dimensions (W x D)

C-MAG MS 4:

C-MAG MS 7:

C-MAG MS 10:

100 x 100mm

180 x 180mm

260 x 260mm

Dimensions (W x D x H)

C-Mag MS 4:

C-Mag MS 7:

C-Mag MS 10:

150 x 260 x 105 mm

220 x 330 x 105 mm

300 x 415 x 105 mm

Weight

C-Mag MS 4 / MS 7 / MS 10:

3kg / 5kg / 6kg

Protection class acc. to DIN EN 60529: IP 21

Tested to DIN EN IEC 61010-1.



Code	Description	Packaging
LLG09720493	Magnetic stirrer C-MAG MS 4 IKAMAG® without heating, glass ceramics set-up plate	1 pz.
LLG09720494	Magnetic stirrer C-MAG MS 7 IKAMAG® without heating, glass ceramics set-up plate	1 pz.
LLG09720495	Magnetic stirrer C-MAG MS 10 IKAMAG® without heating, glass ceramics set-up plate	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

