

**LLG-pH Meter 5****LLG Labware****LLG-pH Meter 5**

The LLG-pH meter 5 is a very easy to use pH/mV/°C handheld meter. The pH value and temperature can be stored directly on the SD card (Excel file) or transferred directly from the pH meter to a PC via the RS-232 interface.

- Gel filled pH electrode with external temperature sensor with 1m cable
- BNC connector
- Automatic temperature compensation (ATC)
- Large, easy-to-read LCD display
- 3 point calibration
- Suitable for all standard SD cards with 1 to 16 GB

Please order batteries separately (6 x 1,5 V AA)

Optional: power cable, PC connection cable, PC software

Scope of supply: pH meter, pH electrode with temperature sensor, buffers pH 4 and pH 7, 4 GB SD-Memory card, transport case

Specifications

Measuring range pH: 0 ... 14

pH resolution: 0.01

pH accuracy: $\pm 0.02 + 2d$

Ambient conditions: 0 ... 50 °C, max. 85 % r.H.

Temperature range electrode: 0 ¼ 60 °C

Dimensions (W x D x H): 68 x 177 x 45 mm

Weight: 490 g

Power supply: 6 x 1.5 V AA batteries / 9 V mains adapter (optional)

IP code: IP 51

Warranty: 3 years (1 year on the electrode)



Code	Description	Packaging
LLG06263691	LLG-pH Meter 5,	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

