

**Microtube Storage Boxes, PP, 50-/100-Well**

Heathrow Scientific

Microtube Storage Boxes, PP, 50-/100-Well

Autoclavable. Compact storage box has angled front slot for easy access to tubes. Rugged, polypropylene box has durable three-point hinges, snapping clasp for secure closure, moulded grid lines on lid, imprinted coordinates on bottom of base and imprinted and raised coordinates on tube wells. Ideal for sample or case study storage. Holds 1.5ml to 2.0ml microtubes. Resist storage temperatures from -80 to +121 °C.



Code	Description	Packaging
LLG06243162	Microtube Storage Boxes, PP, 50-/100-Well, Colour natural, dim. (l x w x h) 141 x 92 x 56 mm, Array 50	1 pz.
LLG09193985	Microtube Storage Boxes, PP, 50-/100-Well, Colour blue, green, purple, yellow, orange, dim. (l x w x h) 141 x 92 x 56 mm, Array 50	5 pz.
LLG06206710	Microtube Storage Boxes, PP, 50-/100-Well, Colour natural, dim. (l x w x h) 141 x 151 x 57 mm, Array 100	1 pz.
LLG09193981	Microtube Storage Boxes, PP, 50-/100-Well, Colour blue, green, purple, yellow, orange, dim. (l x w x h) 141 x 151 x 57 mm, Array 100	5 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

