

**Steam Distillation Unit S****behr****Steam Distillation Unit S**

- Exemplary safety and reliability
- Robust casing of polyurethane
- Distillation time approx.. 3 mins. per sample
- Detection limit 0.1 mg N
- Recovery rate >99.5%
- Reproducibility $\pm 1\%$
- Steam generating capacity adjustable from 10% - 100%
- Extremely simple, menu-controlled operation of the programmable units via one single operating element (behr one-button operation)
- Programmable reaction time
- Serial interface (RS232)
- Filling level monitoring for can set
- Practical quick clamping device equally suited for the "left-handed" analyst



Code	Description	Packaging
LLG06236005	Steam Distillation Unit S, Type KAS 20, Can set for S 1 and S 2	1 pz.
LLG09843816	Steam Distillation Unit S, Type KAS 30, Set of 3 canisters (20 l each) for S 3, incl. float switch	1 pz.
LLG09843818	Steam Distillation Unit S, Type KAS 40, Set of 4 x 20 l containers for S 4, incl. float switch	1 pz.
LLG09843817	Steam Distillation Unit S, Type S 4, Steam distillation unit, fully automatic, like S 3 with automatic extraction H3BO3 of the sample residues and with 99 programs	1 pz.
LLG09843815	Steam Distillation Unit S, Type S 3, Steam distillation unit, partly automated, like S2 with automatic extraction of the sample residues and with 10 programs	1 pz.
LLG09843814	Steam Distillation Unit S, Type S 1, Steam distillation unit, partly automated with automatic addition of NaOH	1 pz.
LLG09843813	Steam Distillation Unit S, Type S 2, Steam distillation unit, partly automated with automatic addition of NaOH and H2O	1 pz.
LLG06270443	Steam Distillation Unit S, Type TB 2, Titration module for connection with S 5	1 pz.
LLG09843803	Steam Distillation Unit S, Type S 5, Steam distillation with automatic addition of NaOH, H2O and H3BO3 with automatic extraction of the sample residues, with 99 programs and a connection option for external titrator	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*



