



Qualitative filter paper, Grade 595 1/2, folded filters

Qualitative filter paper, Grade 595 1/2, folded filters

Medium fast, 4 - 7 μm .

- time saving, because folded filters have a larger surface area
- filter faster than comparable circular filters
- for medium to fine precipitates and general laboratory work
- Weight: 68 g/m²
- Thickness: 0.15 mm



Code	Description	Packaging
LLG09063303	Qualitative filter paper, Grade 595 1/2, folded filters, Diam. 110 mm	100 pz.
LLG09063304	Qualitative filter paper, Grade 595 1/2, folded filters, Diam. 125 mm	100 pz.
LLG09063305	Qualitative filter paper, Grade 595 1/2, folded filters, Diam. 150 mm	100 pz.
LLG09063306	Qualitative filter paper, Grade 595 1/2, folded filters, Diam. 185 mm	100 pz.
LLG09063307	Qualitative filter paper, Grade 595 1/2, folded filters, Diam. 240 mm	100 pz.
LLG09063308	Qualitative filter paper, Grade 595 1/2, folded filters, Diam. 270 mm	100 pz.
LLG09063309	Qualitative filter paper, Grade 595 1/2, folded filters, Diam. 320 mm	100 pz.
LLG09063310	Qualitative filter paper, Grade 595 1/2, folded filters, Diam. 385 mm	100 pz.
LLG09063301	Qualitative filter paper, Grade 595 1/2, folded filters, Diam. 70 mm	100 pz.
LLG09063302	Qualitative filter paper, Grade 595 1/2, folded filters, Diam. 90 mm	



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

