



## Qualitative filter paper, Grade 598½, folded filters

Whatman

### Qualitative filter paper, Grade 598½, folded filters

#### Medium to fast, 8 to 10µm

-thick filter paper

-for fast filtration of coarse precipitates

-Weight: 140g/m<sup>2</sup>

Thickness: 0.32mm



Code	Description	Packaging
LLG09063604	Qualitative filter paper, Grade 598½, folded filters, Diam. 125 mm	50 pz.
LLG09063606	Qualitative filter paper, Grade 598½, folded filters, Diam. 185 mm	50 pz.
LLG09063607	Qualitative filter paper, Grade 598½, folded filters, Diam. 240 mm	50 pz.
LLG09063611	Qualitative filter paper, Grade 598½, folded filters, Diam. 500 mm	50 pz.
LLG09066598	Qualitative filter paper, Grade 598½, folded filters, Diam. 580 x 580 mm	250 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

