



Inserts for LLG-Short Thread Vials ND9/ND10 and LLG-Snap Ring Vials ND11, wide opening

LLG Labware

Inserts for LLG-Short Thread Vials ND9/ND10 and LLG-Snap Ring Vials ND11, wide opening



Code	Description	Packaging
LLG04001547	Inserts for LLG-Short Thread Vials ND9/ND10 and LLG-Snap Ring Vials ND11, wide opening, Nominal capacity 0.3 ml, Cap size ND9/ND10/ND11, Dimensions (ØxH) 6 x 31 mm, Colour Clear, Bottom shape Conical, Version Tip, 12 mm	100 pz
LLG04008196	Inserts for LLG-Short Thread Vials ND9/ND10 and LLG-Snap Ring Vials ND11, wide opening, Nominal capacity 0.2 ml, Cap size ND9/ND10/ND11, Dimensions (ØxH) 6 x 31 mm, Colour Clear, Bottom shape Flat, Material Glass	100 pz
LLG06093247	Inserts for LLG-Short Thread Vials ND9/ND10 and LLG-Snap Ring Vials ND11, wide opening, Nominal capacity 0.25 ml, Cap size ND9/ND10/ND11, Dimensions (ØxH) 5.7 x 29 mm, Colour Clear, Bottom shape Conical, Version Plastic spring	100 pz
LLG07401744	Inserts for LLG-Short Thread Vials ND9/ND10 and LLG-Snap Ring Vials ND11, wide opening, Nominal capacity 0.25 ml, Cap size ND9/ND10/ND11, Dimensions (ØxH) 6 x 31 mm, Colour Clear, Bottom shape Conical, Version Tip, 15 mm	100 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

