


**3-way stopcocks, T-bore, borosilicate glass 3.3****Lenz****3-way stopcocks, T-bore, borosilicate glass 3.3**

Acc. to DIN 12554 and in addition to DIN standard, T-bore, with solid glass plug or PTFE plug and screwthread retaining device. Material: borosilicate glass 3.3, leakage rate to DIN 12540-1.



Code	Description	Packaging
LLG09116031	3-way stopcocks, T-bore, borosilicate glass 3.3, PTFE Key, Stopcock type 14,5 NS, Bore 1.5 mm	
LLG09116011	3-way stopcocks, T-bore, borosilicate glass 3.3, Solid Glass Key, Stopcock type 14,5 NS, Bore 1.5 mm	1 pz.
LLG09116012	3-way stopcocks, T-bore, borosilicate glass 3.3, Solid Glass Key, Stopcock type 14,5 NS, Bore 2.5 mm	1 pz.
LLG09116013	3-way stopcocks, T-bore, borosilicate glass 3.3, Solid Glass Key, Stopcock type 18,8 NS, Bore 2.5 mm	1 pz.
LLG09116033	3-way stopcocks, T-bore, borosilicate glass 3.3, PTFE Key, Stopcock type 18,8 NS, Bore 2.5 mm	
LLG06223066	3-way stopcocks, T-bore, borosilicate glass 3.3, PTFE Key, Stopcock type 18,8 NS, Bore 4 mm	
LLG09116014	3-way stopcocks, T-bore, borosilicate glass 3.3, Solid Glass Key, Stopcock type 18,8 NS, Bore 4 mm	1 pz.
LLG09116015	3-way stopcocks, T-bore, borosilicate glass 3.3, Solid Glass Key*, Stopcock type 24,0 NS, Bore 4 mm	1 pz.
LLG09116035	3-way stopcocks, T-bore, borosilicate glass 3.3, PTFE Key, Stopcock type 24,0 NS, Bore 4 mm	
LLG09116017	3-way stopcocks, T-bore, borosilicate glass 3.3, Solid Glass Key*, Stopcock type 29,2 NS, Bore 6 mm	1 pz.
LLG09116037	3-way stopcocks, T-bore, borosilicate glass 3.3, PTFE Key, Stopcock type 29,2 NS, Bore 6 mm	



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
 Reagents operates with
 a Certified Quality
 Management System

