

**Hotplates, analogue, CB300 and CB500****Stuart****Hotplates, analogue, CB300 and CB500**

Energy regulator control with glass ceramic top plate and 'hot' warning lamp which illuminates when the plate temperature exceeds 50°C. For 230V 50/60 Hz single phase supplies. With BioCote, silver-based, antimicrobial protection. Maximum temperature as indicated.



Code	Description	Packaging
LLG09645316	Hotplates CB500, dim. (W x D x H) 520 x 360 x 130 mm, Hotplate 500 x 300 mm	1 pz.
LLG09645338	Hotplates CB300, dim. (W x D x H) 320 x 370 x 120 mm, Hotplate 300 x 300 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

