

**LLG-Aspirator Bottles, narrow neck, HDPE, with stopcock****LLG Labware****LLG-Aspirator Bottles, narrow neck, HDPE, with stopcock**

With carrying handle and $\frac{3}{4}$ " stopcock. Thread and opening for stopcock provided. Suitable for food. Very strong and durable quality for extended use over time. With gasket in the screw cap for a very tight seal. Screw cap can be attached to the bottle using the integrated loops. The bottles can be used by temperature from -40 °C to +80 °C.

Please use only in upright position.



Code	Description	Packaging
LLG09139933	Aspirator Bottles, with stopcock, 5 l, Neck Int. diam. 52 mm, Neck Ext. diam. 62 mm, Body diam. 160 mm, Height 350 mm, Handle 1	1 pz
LLG09139934	Aspirator Bottles, with stopcock, 10 l, Neck Int. diam. 52 mm, Neck Ext. diam. 62 mm, Body diam. 210 mm, Height 430 mm, Handle 1	1 pz
LLG09139935	Aspirator Bottles, with stopcock, 25 l, Neck Int. diam. 84 mm, Neck Ext. diam. 97 mm, Body diam. 270 mm, Height 570 mm, Handle 2	1 pz
LLG09139936	Aspirator Bottles, with stopcock, 50 l, Neck Int. diam. 84 mm, Neck Ext. diam. 97 mm, Body diam. 350 mm, Height 690 mm, Handle 2	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

