

**LLG-Aspirator Bottles, narrow neck, HDPE, with stopcock****LLG Labware****LLG-Aspirator Bottles, narrow neck, HDPE, with stopcock**

With carrying handle and $\frac{3}{4}$ " stopcock. Thread and opening for stopcock provided. Suitable for food. Very strong and durable quality for extended use over time. With gasket in the screw cap for a very tight seal. Screw cap can be attached to the bottle using the integrated loops. The bottles can be used by temperature from -40 °C to +80 °C.

Please use only in upright position.



Code	Description	Packaging
LLG09139936	Aspirator Bottles, with stopcock, 50 l, Neck Int. diam. 84 mm, Neck Ext. diam. 97 mm, Body diam. 350 mm, Height 690 mm, Handle 2	1 pz
LLG09139933	Aspirator Bottles, with stopcock, 5 l, Neck Int. diam. 52 mm, Neck Ext. diam. 62 mm, Body diam. 160 mm, Height 350 mm, Handle 1	1 pz
LLG09139934	Aspirator Bottles, with stopcock, 10 l, Neck Int. diam. 52 mm, Neck Ext. diam. 62 mm, Body diam. 210 mm, Height 430 mm, Handle 1	1 pz
LLG09139935	Aspirator Bottles, with stopcock, 25 l, Neck Int. diam. 84 mm, Neck Ext. diam. 97 mm, Body diam. 270 mm, Height 570 mm, Handle 2	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

