



Permanent markers, edding 404/400

Permanent markers, edding 404/400

Permanent marker with very fine round tip. Low-odour, smudge-proof, water-resistant, quick drying, lightfast, refillable (only e-400). Replaceable tips. For labelling on almost all surfaces, including glass, metal and plastics. Alcohol-based ink without toluene / xylene.



Code	Description	Packaging
LLG09050421	Permanent markers, edding 404/400, Type 400, Colour Black, Line width 1.0 mm	1 pz
LLG09050451	Permanent markers, edding 404/400, Type 404, Colour Black, Line width 0.75 mm	1 pz
LLG09050423	Permanent markers, edding 404/400, Type 400, Colour Blue, Line width 1.0 mm	1 pz
LLG09050453	Permanent markers, edding 404/400, Type 404, Colour Blue, Line width 0.75 mm	1 pz
LLG09050427	Permanent markers, edding 404/400, Type 400, Colour Brown, Line width 1.0 mm	
LLG09050424	Permanent markers, edding 404/400, Type 400, Colour Green, Line width 1.0 mm	1 pz
LLG09050454	Permanent markers, edding 404/400, Type 404, Colour Green, Line width 0.75 mm	1 pz
LLG09050426	Permanent markers, edding 404/400, Type 400, Colour Orange, Line width 1.0 mm	
LLG09050429	Permanent markers, edding 404/400, Type 400, Colour Pink, Line width 1.0 mm	
LLG09050422	Permanent markers, edding 404/400, Type 400, Colour Red, Line width 1.0 mm	1 pz
LLG09050452	Permanent markers, edding 404/400, Type 404, Colour Red, Line width 0.75 mm	1 pz
LLG09050428	Permanent markers, edding 404/400, Type 400, Colour Violet, Line width 1.0 mm	
LLG09050425	Permanent markers, edding 404/400, Type 400, Colour Yellow, Line width 1.0 mm	



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*



