



Magnetic stirrer C-MAG HS 4 / C-MAG HS 7 / C-MAG HS 10

Magnetic stirrer C-MAG HS 4 / C-MAG HS 7 / C-MAG HS 10

With heating. Ceramic hotplate which offers excellent chemical resistance.

-powerful motor for stirring quantities of up to 5L, 10L or 15L (H₂O).

-fixed safety circuit of 550°C

-Hot Top indicator - hot surface warning to prevent burns

-exact temperature setting via digital display (LED)

-digital error code display

-elevated control panel for protection against spilled liquids

C-MAG HS 7 and C-MAG HS 10:

Connection acc. to DIN 12878 for attaching a contact thermometer, e.g. ETS-D5; enables precise temperature control.

Specifications

Magnetic stirrer HS4 // HS7 // HS10

Stirring quantity (H₂O): 5 l // 10 l // 15 l

Motor rating input / output: 15 / 1.5W

Speed range: 100 ... 1500rpm

Speed display: scale

Heating function

Temperature display:
digital

Heat output: 250W // 1000W // 1500W

Heating rate (1 l H₂O): 2.5 K/min // 5 K/min // 5 K/min

Heating temperature range: 50 ... 500°C

Heat control accuracy: ± 10 K

Fixed safety circuit: 550°C

Control accuracy with sensor

C-MAG HS 4:

C-MAG HS 7 / HS 10:

-

ETS-D5 / ± 0.5 K

Top plate

Material:

glass ceramic

Dimensions: 100 x 100mm // 180 x 180mm // 260 x 260mm

Protection class according to DIN EN 60529: IP21



Code	Description	Packaging
LLG04008087	Magnetic stirrer C-MAG HS 4 S 2 with heating, glass ceramics heating plate with UK plug	1 pz.
LLG04008090	Magnetic stirrer C-MAG HS 7 S 2 with heating, glass ceramics heating plate with UK plug	1 pz.
LLG04008093	Magnetic stirrer C-MAG HS 10 S 2 with heating, glass ceramics heating plate with UK plug	1 pz.
LLG09720490	Magnetic stirrer C-MAG HS 4 with heating, glass ceramics heating plate	1 pz.
LLG09720491	Magnetic stirrer C-MAG HS 7 with heating, glass ceramics heating plate	1 pezzo
LLG09720492	Magnetic stirrer C-MAG HS 10 with heating, glass ceramics heating plate with EU plug	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

